



NetAtlas Installation Guide

Copyright © 2007 by The Hamner Institutes. All rights reserved.

To install NetAtlas plugin on your computer follow the steps below:

1. Install the most recent version of [Cytoscape](#).

NetAtlas is an add-on module for Cytoscape. As a result, proper installation of Cytoscape is required. During installation of Cytoscape, please make a note of the directory that Cytoscape was installed as this information will be needed later on.

NOTE: If Cytoscape has already been installed, please confirm that it is version 2.4 or above and check what directory it has been installed in.

2. Download and save the most recent NetAtlas distribution file (NetAtlas.jar).

3. Copy the NetAtlas distribution file to the Cytoscape plugins folder.

If you installed Cytoscape to the default directory, the location of the plugin folder will be:

A. Windows Systems: C:\Program Files\Cytoscape_v2.4.x\plugins

B. Macintosh System: Macintosh HD\Applications\Cytoscape_v2.4.x\plugins

C. Linux System: /usr/local/Cytoscape_v2.4.1/plugins

4. Restart Cytoscape

The NetAtlas plugin will be available under the 'Plugins' pulldown menu.